

Cryotainer B.V.
Fokkerstraat 507
3125 BD Schiedam
Phone: +31 (0)10 262 30 53

P.O. Box 130
3100 AC Schiedam
www.cryotainer.com



Datasheet FT 1700 (Flying Tainer)



Dimensions l x h x b	mm	12.200 x 3.650 x 2.500
Maximum working pressure	bar	10,0 or 22,0 (depending on storage container)
Maximum flow	Nm ³ /h	1.700 (2 x 850) for 8 hours per day maximum
Maximum outlet temperature nitrogen	°C	15 ⁰ till 20 ⁰ below ambient temperature
Total weight	kg	38.000
Suitable for cryogenic gasses		LIN; LOX*
Storage capacity	ltr	19.700
Outlet connection		Both sides 2" weco fig. 100

* For use of LOX container and vaporizer should be purged before use